


# Advanced Hemodynamic Monitoring with the Edwards Swan-Ganz Catheter

The advanced technology Swan-Ganz catheter, using the same functionality as the standard Swan-Ganz catheter, provides the ability to continuously monitor the balance between oxygen delivery and consumption and investigate the root cause of an imbalance through analysis of the components of stroke volume (preload, afterload, and contractility).




## ADVANTAGES OF THE ADVANCED TECHNOLOGY SWAN-GANZ CATHETER AS COMPARED TO THE STANDARD SWAN-GANZ CATHETER

- Maximum amount of diagnostic information with same invasive procedure
- Continuous assessment of  $DO_2 / VO_2$  balance with  $SvO_2$  monitoring
- Continuous assessment of adequacy of CO by assessing  $DO_2 / VO_2$  balance with  $SvO_2$  monitoring
- Continuous assessment of components of stroke volume (preload, afterload, and contractility) (RVEDI, SVR, RVEF and SVI)
- Increased accuracy of cardiac output calculations, elimination of ventilator cycle and thermal noise effect
- Mitigation of:
  - User error in association with wedge procedure/calculation through automated alternative preload parameter (RVEDI)
  - Inappropriate preload assessment secondary to changes in ventricular compliance affecting PAD or PAOP
  - Pulmonary artery rupture potential associated with wedge procedure by providing automated preload parameter (RVEDI)
  - Inappropriate therapy due to miscalculation of PAOP by using automated preload parameter (RVEDI)
  - Iatrogenic infection risk from bolus injections
  - Cardiac output error with CCO automation through elimination of bolus cardiac output user error



### TROUBLESHOOTING CCO AND CCOMBO CATHETERS

Problem/Displayed Message	Possible Cause	Suggested Action
<b>CCO Measurements Do Not Agree With Bolus CO Measurements</b>	<ul style="list-style-type: none"> <li>• Catheter position</li> <li>• Instrument configuration</li> <li>• Faulty thermistor or injectate probe connection</li> <li>• Unstable baseline temperature affecting bolus CO measurements</li> </ul>	<ul style="list-style-type: none"> <li>• Verify proper catheter position:                             <ul style="list-style-type: none"> <li>- Balloon inflation 1.5 cc</li> <li>- Check for PAOP tracing</li> <li>- Transduce injectate (RA/blue) lumen</li> </ul> </li> <li>• Verify that injectate volume and catheter size OR computation constant have been correctly selected</li> <li>• Check injectate probe and catheter connection</li> <li>• Assess factors affecting PA temperature stability:                             <ul style="list-style-type: none"> <li>- Rapid bolus infusions</li> <li>- 60 seconds between injections</li> <li>- Shivering</li> <li>- Patient movement</li> <li>- Fighting ventilator</li> </ul> </li> </ul>
<b>*Check Thermal Filament Position*</b>	<ul style="list-style-type: none"> <li>• Flow around thermal filament may be reduced</li> </ul>	<ul style="list-style-type: none"> <li>• Verify proper catheter position</li> <li>• Verify free-floating state of catheter</li> </ul>
<b>*Warm Injectate* or *Injectate too Warm Check Probe*</b>	<ul style="list-style-type: none"> <li>• Injectate temperature within 8°C of blood temperature</li> <li>• Injectate temperature &gt;30°C</li> </ul>	<ul style="list-style-type: none"> <li>• Use cooler injectate fluid:                             <ul style="list-style-type: none"> <li>- Maintain ice slush solution</li> <li>- If pre-filled syringes, remove from cold source and use within 15 seconds</li> </ul> </li> <li>• Check injectate probe connection</li> <li>• Replace injectate temperature probe</li> </ul>
<b>Loss of Power to Instrument</b>	<ul style="list-style-type: none"> <li>• Power surge</li> <li>• Power outage</li> </ul>	<ul style="list-style-type: none"> <li>• Press <math>SvO_2</math></li> <li>• Press TRANSPORT</li> <li>• Press RECALL (Note: Calibration data must be less than 24 hours old)</li> </ul>
<b>*SQI = 4*</b>	<ul style="list-style-type: none"> <li>• Catheter too distal</li> <li>• Low blood flow at catheter tip or catheter tip against vessel wall</li> <li>• Change in Hgb/Hct values</li> <li>• Catheter kinked or damaged</li> </ul>	<ul style="list-style-type: none"> <li>• Verify proper catheter position</li> <li>• Verify patency of catheter</li> <li>• Update Hgb/Hct values using UPDATE function</li> <li>• Check catheter for kinking and recalibrate</li> <li>• Replace catheter if values required and recalibrate</li> <li>• After all above troubleshooting measures are utilized, press OPTICAL RESET</li> </ul>
<b>*Red/IR Transmit*</b>	<ul style="list-style-type: none"> <li>• Optical module contamination/Optical module damage</li> </ul>	<ul style="list-style-type: none"> <li>• Disconnect catheter from OM</li> <li>• Clean optical connector on catheter and OM</li> <li>• Reconnect. If error message disappears, do INVIVO calibration</li> <li>• If error message remains, try different optical module</li> <li>• If error message remains, change catheter</li> </ul>



### RVEDI / RVEF

1. Parameters Attained with Vigilance II Monitor

- End-Diastolic Volume (EDV): The volume of blood in the ventricle at the end of the diastole.  $EDV = SV/EF$   
Normal RVEDI: 100-160 mL  
Normal RVEDVI: 60-100 mL/m<sup>2</sup>
- Ejection Fraction (EF): The percentage of blood ejected from the ventricle each beat  
 $EF = \frac{EDV - ESV}{EDV}$  or  $\frac{SV}{EDV}$   
Normal RVEF: 40 - 60%
- Stroke Volume (SV): The volume of blood ejected from the ventricle in each beat  
 $SV = CO / HR \times 1000$   
Normal SV: 60 - 100 mL  
Normal SVI: 33 - 47 mL/beat/m<sup>2</sup>

(Note: As with all measurements in hemodynamic monitoring, the absolute number is not as important as trends and changes in response to therapy)



2. Goal of RV Volumetric Measurements

- Optimize RV efficiency
- Optimize the relationship between EDV and SV
- A. In an efficient state, an increase in preload (EDV) will result in an increase in stroke volume (SV)
- B. Prior to reaching the flat part of the curve, an increase in preload (EDV) will increase SV while not causing a decrease in ejection fraction
- C. On the flat part of the curve, a further increase in preload (EDV) will not result in an increase in SV

At this point, a further increase in volume may:

- Decrease oxygen supply
- Increase oxygen demand
- Decrease left ventricular compliance

Therapy should be directed at increasing contractility or reducing afterload




### APPLICATIONS AND CONTRAINDICATIONS

**Clinical applications for Swan-Ganz pulmonary artery catheters:**

- Intra-abdominal hypertension
- Patients at risk for acute right ventricular dysfunction
- ARDS
- Extensive burns
- Cardiac surgery
- Significant cardiac tamponade
- Significant cardiomyopathy
- Significant constrictive pericarditis
- Drug intoxication
- Severe eclampsia
- Significant intra- or extra-vascular fluid shifts
- At risk for hemorrhage
- Intra- and post-operative high-risk surgery management
- Patient on intra-aortic balloon counterpulsation
- Complex liver resections
- Liver transplantation
- Complex lung resection
- Complex myocardial infarctions
- Pulmonary edema
- Pulmonary embolism
- Pulmonary hypertension
- Acute renal failure
- Severe sepsis
- Presence of or at risk for cardiogenic shock
- Presence of or at risk for distributive shock
- Presence of or at risk for hemorrhagic shock
- Shock of unknown etiology
- Shock unresponsive to attempts at resuscitation
- Severe trauma
- Ventilator effects on hemodynamics

**Relative contraindications for Swan-Ganz pulmonary artery catheterization:**  
(There are no absolute contraindications to the use of a pulmonary artery catheter; risk/benefit must be assessed for each patient)


- Left bundle branch block
- Patients with tricuspid or pulmonic heart valve replacements
- Presence of endocardial pacing leads
- Lack of appropriate clinical skills or infrastructure to insert and/or support the use of a pulmonary artery catheter
- Heparin coated catheters in patients with known sensitivity to heparin



### CATHETER INSERTION DISTANCE MARKINGS

Location	Distance to VC/RA Junction	Distance to PA
Internal Jugular	15-20 cm	40-55 cm
Subclavian Vein	10-15 cm	35-50 cm
Femoral Vein	30 cm	60 cm
Right Antecubital Fossa	40 cm	70 cm
Left Antecubital Fossa	50 cm	80 cm

Note: Catheter markings occur every 10 cm and are denoted by a thin black line. 50 cm markings are denoted by a thick black line. Catheter must exit introducer sheath before inflating balloon.



**References:**

1. Alspach JG (ed.). Core Curriculum for Critical Care Nursing. 6th ed. St. Louis: Saunders Elsevier; 2006.
2. Daily EK, Schroeder JS. Techniques in bedside hemodynamic monitoring. 5th ed. St. Louis: Mosby; 1994.
3. Darovic, GO. Hemodynamic monitoring: invasive and noninvasive clinical application. 3rd ed. Philadelphia: Saunders; 2002.
4. Headley, JM. Invasive Hemodynamic Monitoring: Applying Advanced Technologies. Critical Care Nursing Quarterly. 1998;21:3:73-84.
5. Headley, JM. Puzzled by Continuous Cardiac Output Monitoring? Nursing '97. 1997;32aa - 32dd.
6. Headley, JM. Special Pulmonary Procedures; Continuous Mixed Venous Oxygen Saturation Monitoring. In: Lynn-McHale DJ, Carlson KK. AACN Procedure Manual for Critical Care. 5th Edition. Philadelphia.: W.B. Saunders Company; 2005.
7. Imperial-Perez F, McRae M, Gawlinski A, Keckeisen M, Jesurum J. AACN protocols for practice: Hemodynamic Monitoring. 1998.
8. Leeper B. Monitoring right ventricular volumes: a paradigm shift. AACN Clin Issues. 2003 May; (14):208-19.
9. Mims BC, Toto KH, Luecke LE, Roberts MK, Brock JD, Tyner TE. Critical care skills: a clinical handbook. 2nd ed. St. Louis: Saunders; 2004.
10. Perret C, Tagan D, Feihl F, Marini JJ. The pulmonary artery catheter in critical care. Cambridge: Blackwell Science Inc.; 1996.
11. Pinsky MR, Vincent JL. Let us use the pulmonary artery catheter correctly and only when we need it. Crit Care Med. 2005 May;33(5):1119-22.
12. Thelan LA, Davie JK, Urden LD, Lough ME. Critical care nursing: diagnosis and management. 2nd ed. St. Louis: Mosby; 1994.
13. Vincent JL, Pinsky MR, Sprung CL, Levy M, Marini JJ, Payen D, Rhodes A, Takala J. The pulmonary artery catheter: in medio virtus. Crit Care Med. 2008 Nov;36(11):3093-6.
14. Woods SL, Froelicher ESS, Motzer SU, Bridges EJ. Cardiac nursing. 5th ed. Philadelphia: Lippincott Williams & Wilkins; 2005.
15. Zink W, Nöhl J, Rauch H, Bauer H, Desimone R, Martin E, Böttiger BW. Continuous assessment of right ventricular ejection fraction: new pulmonary artery catheter versus transoesophageal echocardiography. Anaesthesia. 2004 Nov;59(11):1126-32.

William McGee is a paid consultant of Edwards Lifesciences.

Rx only. See instructions for use for full prescribing information.

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

Edwards, CCombo, and Vigilance II are trademarks of Edwards Lifesciences Corporation, Edwards Lifesciences, the stylized E logo, CCombo V, Swan-Ganz, Vigilance and VIP are trademarks of Edwards Lifesciences Corporation and are registered in the United States Patent and Trademark Office.

© 2009 Edwards Lifesciences LLC. All rights reserved. AR04352

